

Link to the product: <https://sforne.com/moressa-osian-set-premium-turquoise-silicone-p-4703.html>

## MORESSA - OSIAN SET PREMIUM TURQUOISE SILICONE



Price	<b>50.06 £</b>
Number	<b>70137545</b>
Producer code	<b>FETISHID-213249</b>
EAN	<b>8425402156018</b>

### Product description

Moressa OSIAN Chinese balls are covered in 100% high quality silicone. They have an elegant design in three colors as well as an incomparable soft and very comfortable velvet touch.

OSIAN SET is the ideal pelvic floor training complement, especially for beginner women and/or with womb inclinations. It is inserted like a tampon and, like it, OSIAN ONE is almost imperceptible. During everyday movements, the constant and silent internal rotations generate involuntary muscle contractions that strengthen the pelvic floor muscles.

OSIAN SET helps to tone the pelvic floor muscles responsible for sexual pleasure, they have a perfectly dosed weight that guarantees a sensational pelvic floor training session, to reduce incontinence, intensify orgasm and consequently improve sexual health.

For greater effectiveness, short and regular sessions are recommended. Start with OSIAN ONE for a maximum of 15 minutes a day when you are on the move, once you can retain them easily continue with OSIAN TWO, and repeat the process, when you manage to retain them continue with OSIAN THREE repeating the last step again and congratulations! You have achieved perfect pelvic floor strengthening with Amoressa Set.

Moressa balls are completely smooth, making insertion easier and producing a feeling of comfort. Easy to use and clean, it is a perfect choice for impeccable pelvic muscle health.

- 100% hypoallergenic silicone structure.
- Weight; Osian One; 51 Gr / Osian Two; 101 Gr / Osian Three 85 gr
- Measures; Osian One / 34mm – Osian Two / 35mm – Osian Three 34 mm
- Once you have them inside you, your pelvic muscles will start working trying to keep the balls in place.
- The rotation of the ball inside stimulates and strengthens the pelvic floor muscles during the movements of daily life.
- Increases sensitivity during sex, both for WOMEN and MEN!
- Midwives and gynecologists recommend it as a method of recovering the pelvic floor after childbirth.
- Increases the results of postpartum gymnastics
- Prevents urinary incontinence and uterine abruptions
- Faint

### Use tips

Moressa recommends the use of water-based lubricants for easy use and insertion of all silicone products. Do not use silicone-based lubricants, oils or creams with silicone products as they can degrade the material, limiting its operation and useful life.

### Cleaning

Thoroughly clean the toy before and after each use with warm water and neutral soap.

### Storage

---

Store Moressa Osian balls in a dust-free place and keep it out of contact with other toys made of different materials (TPR, PVC jelly, etc.) Avoid exposing it to sunlight and extremely high temperatures.

- 2 years warranty

## Kegel Guide

### **Benefits of Kegel Training**

Stronger pelvic floor muscles improve bladder control, enhance sexual sensation (for both partners), and speed postpartum recovery.

Most users notice improvements within 4-6 weeks of regular (daily) use.