Link to the product: https://sforne.com/moressa-osian-three-premium-lilac-silicone-p-4689.html



MORESSA - OSIAN THREE PREMIUM LILAC SILICONE

Price	17.51 €
Number	70137541
Producer code	FETISHID-213245
EAN	8425402155981

Product description

Moressa OSIAN Chinese balls are covered in 100% high quality silicone. They have an elegant design in three colors as well as an incomparable soft and very comfortable velvet touch.

OSIAN THREE is the ideal pelvic floor training supplement, especially for beginner women and/or with womb inclinations. It is inserted like a tampon and, like it, OSIAN ONE is almost imperceptible. During everyday movements, the constant and silent internal rotations generate involuntary muscle contractions that strengthen the pelvic floor muscles.

OSIAN THREE helps to tone the pelvic floor muscles responsible for sexual pleasure, they have a perfectly dosed weight that guarantees a sensational pelvic floor training session, to reduce incontinence, intensify orgasm and consequently improve sexual health.

For greater effectiveness, short and regular sessions are recommended. Once the pelvic floor muscles have been strengthened with OSIAN THREE, repeat the procedure with OSIAN ONE, TWO and THREE periodically every 6 months.

Moressa balls are completely smooth, making insertion easier and producing a feeling of comfort. Easy to use and clean, it is a perfect choice for impeccable pelvic muscle health.

- 100% hypoallergenic silicone structure.
- Weight; 85g
- Measures; 34mm diameter
- Once you have them inside you, your pelvic muscles will start working trying to keep the balls in place.
- The rotation of the ball inside stimulates and strengthens the pelvic floor muscles during the movements of daily life.
- Increases sensitivity during sex, both for WOMEN and MEN!
- Midwives and gynecologists recommend it as a method of recovering the pelvic floor after childbirth.
- Increases the results of postpartum gymnastics
- Prevents urinary incontinence and uterine abruptions
- Faint

Use tips

Moressa recommends the use of water-based lubricants for easy use and insertion of all silicone products. Do not use silicone-based lubricants, oils or creams with silicone products as they can degrade the material, limiting its operation and useful life.

Cleaning

Thoroughly clean the toy before and after each use with warm water and neutral soap.

Storage

Store Moressa Osian balls in a dust-free place and keep it out of contact with other toys made of different materials (TPR, PVC

jelly, etc.) Avoid exposing it to sunlight and extremely high temperatures.	
• 2 years warranty	
, ,	
generated by the shortCold	