

Link to the product: <https://sforne.com/packaging-available-in-es/en/fr/de/pt/it/ro/el/ru/hu/>



PACKAGING AVAILABLE IN: /es/en/fr/de/pt/it/ro/el/ru/hu/

| | |
|---------------|------------------------|
| Price | 12.95 Euro |
| Number | 70140008 |
| Producer code | FETISHID-229297 |
| EAN | 8435565930837 |

Product description

Our BIOAQUA lubricant is the perfect choice for natural pleasure: no spin, just frictionless sensations with a light touch. It can be used during sex, but thanks to its natural water-based formulation, it is also gentle with sex toys.

What can BIOAQUA lubricant provide?

The BIOAQUA version of the water-based lubricant, thanks to its soft and 100% natural formula, offers a silky skin sensation that makes the stimulation especially intense. It is suitable for both women and men or couples and always provides greater sliding pleasure. The formula without oil or greasy components is also compatible with condoms.

This water-based lubricant is 100% natural, safe for toys and pleasant for the body!

BIOAQUA is also great for use with condoms or just skin-to-skin. Unlike many other water-based lubricants, BIOAQUA does not leave a residue. Also, its transparent bottle looks great on your nightstand.

Please note: The use of non-water based lubricants, oils and creams can cause irreversible damage to the control unit and the entire toy.

Instructions for use: apply the amount of lubricant you want on the desired parts of the body and repeat it as many times as you want.

BIOAQUA lubricants are compatible with Latex condoms.

- Suitable for oral sex
- Ideal for use with sex toys
- 100% European.

Please note: Using quality non-water based lubricants, oils and creams can cause irreversible damage to the control unit and the entire toy.

- 100% natural ingredients.
- Dermatologically tested
- High performance
- Long sliding capacity
- Suitable for use with latex condoms
- Water soluble
- Neutral to taste
- Colorless and odorless
- Does not contain oil or perfumes
- Approved as a medicinal product with the CE mark

