

Link to the product: <https://sforne.com/satisfyer-woman-disinfectant-spray-300-ml-p-7201.html>



SATISFYER - WOMAN DISINFECTANT SPRAY 300 ML

| | |
|---------------|------------------------|
| Price | 16.22 £ |
| Number | 70140073 |
| Producer code | FETISHID-231297 |
| EAN | 4019514306603 |

Product description

Product information "Satisfyer Woman Disinfectant Spray"

- This alcohol-free disinfectant spray provides simple and reliable protection against bacteria and viruses. Thanks to its convenient spray head, cleaning your pleasure products are easy, gentle on the material, and effective - all at the same time.

Why should I use the disinfectant spray for women?

- Your Satisfyer is designed to give you intense and breathtaking orgasms on a regular basis. To ensure that the delicate, super-soft silicone retains its high quality for a long time, it's important to always clean and maintain the surface thoroughly. The Satisfyer disinfectant spray for women is specially formulated for your Satisfyer and carefully removes a wide variety of fungi, bacteria, and viruses. The alcohol-free Satisfyer Woman Disinfectant Spray is extra strong but still gentle on the material. Do not use a conventional disinfectant when cleaning, as this can attack the soft surface of your Satisfyer, make it porous, or damage your skin if you use your device again immediately after cleaning. Our Satisfyer sexual wellness device cleaners are harmless to your health, but should always be used directly on the product itself and not on your skin.

How do I use the disinfectant spray for women?

- Thanks to the practical spray attachment you can easily wet your Satisfyer all around with disinfectant. First clean your pleasure-bringer with lukewarm water and soap and dry it with a dust-free cloth. Thanks to the waterproof (IPX7) finish of Satisfyer products, you don't have to worry when cleaning it with water. Then spray the product generously with the disinfectant and let it take effect for 1-2 minutes. Just wipe it dry with a lint-free cloth and your Satisfyer is ready for your next orgasm.

Properties "Satisfyer Woman Disinfectant Spray"

- Content: 300 ml
- Weight: 340 g
- Contents: 100g of solution contains 0,24g Didecyldiemethylammoniumchlorid, 0,24g Alkyldimethylbenzylammoniumchlorid, 0,24g Alkyldimethylethylammoniumchlorid

Cleaning Guide

Cleaning Guide by Material

Silicone: soap & water, boiling, or bleach solution. Glass/Metal: soap & water, boiling, or dishwasher (top rack). TPE/TPR: soap & water only (no boiling). ABS plastic: soap & water only.

Never use silicone lube on silicone toys - use water-based lube instead.