

Link to the product: <https://sforne.com/satisfyer-yoni-power-1-training-set-white-p-5214.html>



SATISFYER - YONI POWER 1 TRAINING SET WHITE

Price	24.79 €
Number	70138066
Producer code	FETISHID-229932
EAN	4049369043835

Product description

You can use the Satisfyer Yoni Power 1 to train your pelvic floor muscles in just 15 minutes a day to experience even more intense orgasms. The different sizes of these Ben Wa balls allow a slow, gradual increase.

Product information "Yoni Power 1"

- Pelvic floor training strengthens the muscles for more intense orgasms
- Made from super soft, skin-friendly medical-grade silicone that is smooth to the touch and extremely hygienic
- A retaining ring enables safe use
- 15-year guarantee
- Training set to strengthen the muscles of the pelvic floor
- More fun for both men and women during sex thanks to muscular constrictions
- Safe use thanks to the return strap
- Medical-grade silicone
- Easy to clean

Satisfyer Yoni Power 1: Pelvic floor training with Ben Wa balls

With the Satisfyer Yoni Power 1 you can strengthen your pelvic floor and enjoy even more intense orgasms. In contrast to pelvic floor exercises with other Ben Wa balls, the focus of the Yoni Power 1 is not just on the weight, but above all on the size of the balls: You start with the largest and heaviest ball (3 cm diameter, 74 g) and work your way through the medium (2.5 cm diameter, 46 g) to the smallest and lightest ball (2 cm diameter, 22 g). The more frequently you do the pelvic floor training, the tighter your muscles become, so you can eventually hold the ball with the smallest circumference. You only have to use the balls for 15 minutes a day to achieve an effect that not only you, but also your sweetheart will notice during vaginal intercourse. The ribbed structure of the round balls also stimulates you during training.

Other great things about the Satisfyer Yoni Power 1 Ben Wa balls

The inner balls of the Satisfyer Yoni Power 1 are made of stainless steel and are encased in soft, medical-grade silicone, which makes pelvic floor training particularly gentle and hygienic. Thanks to their rounded shape and seamless processing, the Ben Wa balls are particularly easy to insert and just as easy to remove thanks to the silicone loops – making the Satisfyer Yoni Power 1 ideal for beginners. The Satisfyer Ben Wa balls can also be used for pelvic floor training after giving birth or to prevent bladder weakness. Use the Ben Wa balls in conjunction with your favorite air-pulse vibrator to further increase your sensations in sexy solo moments. After use, you can clean them with warm water and mild soap and disinfect them with a sexual wellness device cleaner.